

SF Surgical Gloves Latex



Product Description

SF Surgical Gloves Latex for single use, sterile and non-sterile, powdered and powder-free, non-toxic, non-pyrogenic.

Quality

Quality management system according to the requirements of DIN EN ISO 13485

CE Certification: Medical device class 2a according to the requirements of MDD 93/42/EEC

Standart conformity

EN 455 1-4, ISO 10282

Product features


- for surgical operations and examinations
- textured or smooth surface
- anatomically shaped (curved fingers)
- with cuff and cuff edge
- made of natural latex (natural rubber latex)
- powdered or powder-free
- individually and sterile blistered or sealed, in intermediate packaging
- Shelf life 5 years
- Gamma sterilization

Physical parameters

Size	Length, in mm, min	Width, in mm, min	Thickness single wall, in mm, min	Force at break	Elongation
5,0	250	67±4	all sizes finger / palm min 0,10** (smooth surface) min 0,13** (textured surface)	<i>before Aging</i>	
5,5	250	72±4		min 9,0 N* / min 12,5 N**	min 700%**
6,0	260	77±5			
6,5	260	83±5		<i>after Aging</i>	
7,0	270	89±5		min 9,0 N* / min 9,5 N**	min 550%**
7,5	270	95±5			
8,0	270	102±6			
8,5	280	108±6			
9,0	280	114±6			
9,5	280	121±6			

*acc. to EN 455-2, **acc. to ISO 10282

Article description	Size	Color	Number of pieces per package	Article number
Surgical Gloves Latex sterile powdered	5,0	natural white	50 pairs	22.100.000
Surgical Gloves Latex sterile powdered	5,5	natural white	50 pairs	22.100.001
Surgical Gloves Latex sterile powdered	6,0	natural white	50 pairs	22.100.002
Surgical Gloves Latex sterile powdered	6,5	natural white	50 pairs	22.100.003
Surgical Gloves Latex sterile powdered	7,0	natural white	50 pairs	22.100.004
Surgical Gloves Latex sterile powdered	7,5	natural white	50 pairs	22.100.005
Surgical Gloves Latex sterile powdered	8,0	natural white	50 pairs	22.100.006
Surgical Gloves Latex sterile powdered	8,5	natural white	50 pairs	22.100.007
Surgical Gloves Latex sterile powdered	9,0	natural white	50 pairs	22.100.008
Surgical Gloves Latex sterile powdered	9,5	natural white	50 pairs	22.100.009
Surgical Gloves Latex sterile powder-free	5,0	natural white	50 pairs	22.101.000
Surgical Gloves Latex sterile powder-free	5,5	natural white	50 pairs	22.101.001
Surgical Gloves Latex sterile powder-free	6,0	natural white	50 pairs	22.101.002
Surgical Gloves Latex sterile powder-free	6,5	natural white	50 pairs	22.101.003
Surgical Gloves Latex sterile powder-free	7,0	natural white	50 pairs	22.101.004
Surgical Gloves Latex sterile powder-free	7,5	natural white	50 pairs	22.101.005
Surgical Gloves Latex sterile powder-free	8,0	natural white	50 pairs	22.101.006
Surgical Gloves Latex sterile powder-free	8,5	natural white	50 pairs	22.101.007
Surgical Gloves Latex sterile powder-free	9,0	natural white	50 pairs	22.101.008
Surgical Gloves Latex sterile powder-free	9,5	natural white	50 pairs	22.101.009
Surgical Gloves Latex Anatomical non-sterile powdered	5,0	natural white	50 pairs	22.102.000
Surgical Gloves Latex Anatomical non-sterile powdered	5,5	natural white	50 pairs	22.102.001
Surgical Gloves Latex Anatomical non-sterile powdered	6,0	natural white	50 pairs	22.102.002
Surgical Gloves Latex Anatomical non-sterile powdered	6,5	natural white	50 pairs	22.102.003
Surgical Gloves Latex Anatomical non-sterile powdered	7,0	natural white	50 pairs	22.102.004
Surgical Gloves Latex Anatomical non-sterile powdered	7,5	natural white	50 pairs	22.102.005
Surgical Gloves Latex Anatomical non-sterile powdered	8,0	natural white	50 pairs	22.102.006
Surgical Gloves Latex Anatomical non-sterile powdered	8,5	natural white	50 pairs	22.102.007
Surgical Gloves Latex Anatomical non-sterile powdered	9,0	natural white	50 pairs	22.102.008
Surgical Gloves Latex Anatomical non-sterile powdered	9,5	natural white	50 pairs	22.102.009

 Product markings:

