

SF Venox Peripheral IV Catheter



Product Description

Peripheral IV Catheter for continuous peripheral intravenous infusion, for single use, sterile, non-toxic, pyrogen free

Quality

- Quality Management System according to DIN EN ISO 13485
- CE marking: Medical Device Class 2a according to MDD 93/42/EEC

Standard conformity:

- DIN EN ISO 10555-1, DIN EN ISO 10555-5

Product features

- Special tapered and kink resistant radiopaque catheter for easy insertion with optimal flow rate
- Special Flexible wings to provide easy and secure fixation
- Manufactured from tested bio-compatible materials offering longer indwelling time
- Rounded grip on the needle hub for better support during insertion
- Injection port with non-return silicon valve for intermittent medication
- Transparent flashback chamber for easy visual medication
- Standard colour coding according to catheter size

Sizes	14G, 16G, 17G, 18G, 20G, 22G, 24G, 26G
Material	PUR or FEP
Configuration	With flashback chamber, with or without safety features, with or without port, with small or suturable or without wings, with or without snap fit cap, with or without three-way-stop-cock, with or without dual flashback, with or without hydrofobic bacteria retention filter
Packaging	Paper Blister or Tyvek Blister
Shelf life	5 years
Sterilization	Ethylene oxide

Articles**SF Venox Peripheral IV Catheter in different modifications****With safety features, material: PUR**

- SF Venox PUR Safety Peripheral IV catheter with port, snap fit cap and wings
- SF Venox PUR Safety Peripheral IV catheter with integrated 3-way stop cock
- SF Venox PUR Safety Peripheral IV Catheter with port, snap fit cap, wings and dual flashback

Without safety features, material: PUR

- SF Venox PUR Peripheral IV Catheter with port, snap fit cap and wings
- SF Venox PUR Peripheral IV Catheter with integrated 3-way-stop-cock
- SF Venox PUR Peripheral IV Catheter without port and wing
- SF Venox PUR Peripheral IV Catheter with port and wings
- SF Venox PUR Peripheral IV Catheter without port and with wings
- SF Venox PUR Peripheral IV Catheter without port and with small wings
- SF Venox PUR Peripheral IV Catheter with port, snap fit cap, wings and dual flashback

Without safety features, material: FEP

- SF Venox Peripheral IV Catheter with port, snap fit cap and wings
- SF Venox Peripheral IV Catheter with integrated 3-way-stop-cock
- SF Venox Peripheral IV Catheter without port and wing
- SF Venox Peripheral IV Catheter with port and wings
- SF Venox Peripheral IV Catheter without port and with wings

More information about the different variants you are interested in can be obtained on request

Product markings

STERILE | EO